



We listened to infusion nurses

Short staffing



reduced set-up time

Shift handoffs



easy labeling system

Reduced
complexity



multi-fluid management
without electronic devices

Optimizing KVO



AutoStart™ flushes automatically
US Pat. 9,352,080

Flow selector
Pat. pending



ready flush of syringe/mini IV bag
eliminates dead volume loss

Hazardous
protection



fully closed system protects
you and your patient



Redefining Gravity-Fed IV Infusions

Delivering safety, simplicity, and labor savings
while assuring vein patency

Safer Infusion System's enhanced AutoStart™ burette functionality delivering IV infusion excellence for every infusion environment

Our mission: Providing IV infusion systems that optimize caregiver ability to deliver all “10 rights” of medication administration

- 1) Right Patient
- 2) Right Medication
- 3) Right Dose
- 4) Right Route
- 5) Right Time & Frequency
- 6) Right Education
- 7) Right Reason
- 8) Right to Refuse
- 9) Right Response
- 10) Right Documentation



Initially brought to market as the “AutoStart Burette” in Australia, nursing feedback triggered an extensive redesign receiving additional patent approval delivering enhanced functionality, simplicity, and safety in a low-cost consumable.



Dr. Michael Monsour, who received the Order of Australia medal, is Chairman of Analytica and now Safer Infusion Systems, Inc. In keeping with his exemplary medical services, he was the driving force behind this IV infusion technology solving two problems long associated with traditional burettes.

“First, the drip does not have to be turned off when injecting medicine as with traditional burettes. This eliminates the possibility of the drip not being re-established after the drug dose. Additionally, the caregiver does not need to return to the bedside to determine if all of the medicine has been delivered to re-establish the drip. AutoStart’s shuttle valve “restarts” the drip once the medicine has passed through the burette. This saves time for busy caregivers,” advises Dr. Monsour.

“Second, the shuttle valve float action prevents air from leaving the burette eliminating any air embolism event. After flushing the medicine dose and exhaustion of a primary IV bag, the float’s weight seals the exit nozzle. This tight seal retains a small fluid volume in the burette base.”

Relying on a gravity feed eliminates complexity, training, capital costs, and maintenance with infusion pumps making the SIS ideal for hospitals and infusion centers while paving a smoother, safer, and lower-cost path for hospital-to-home IV infusions. The many capabilities of the single SKU Safer Infusion System will shine for tight spaces and mobility demanded for military deployment and emergency transport vehicles.



US 9,352,080 B2

Time savings

Hang one large Primary IV bag (typ. 1,000ml) to prime SIS system + support AutoStart slow and steady catheter flow delivery

Color-coded Primary & Secondary IV bag spikes




Safety enhancement

Consistency of gravity

SIS designed for **18" to 24"** of static **head height** (H_{TH}) above the patient's **head IV site** with a standard IV pole

$$H_{TH} = H_{BL} + H_{DCL}$$





Simpler IV rate control

Focus on drip rate control by roller clamp below the drip chamber - positioned in the middle of the Caregiver Visual Operating Window



AutoStart flows down the central tube in the burette automatically controlled by float shuttle valve

Maintain vein patency with uninterrupted IV flow with minimal caregiver intervention



Protection against CRBSI and occlusion from fibrin build-up



Ready for multiple IV bags




Ready Drug Port syringe flushing


Ready to provide

Micro or Macro Drip protocol

Gentle flow gravity drip best maintains vein patency

No pump-related micro-pistoning of catheter in the vein





Set Up of Primary IV Bag

Connect an IV Bag. Fluid will flow into and fill the chamber and stop at approximately 10 milliliters

The White Roller Clamp / Valve below the drip chamber controls the drip rate for infusion flow to the patient's catheter

Medication Delivery

Secondary IV Bag - Dosed

Turn **Flow Selection Knob** to the North

Open **Pink Roller Clamp** – Dispense contents of dosed **Secondary IV Bag** into Chamber

Medication flow rate is managed by the White Roller Clamp establishing the prescribed drip rate in the **Drip Chamber**

Upon exhaustion of the **Secondary IV Bag** fluid in the chamber, **AutoStart** mode commences to steadily flush the IV line supporting vein patency at the catheter site

Secondary IV Bag – Dosing of Undosed Bag

Turn **Flow Selection Knob** to the North

Open **Pink Roller Clamp** – to prime fluid down to the Chamber & turn **Flow Selection Knob** to the East

Attach injector for medicinal load to **Needless Drug Port Connector** & push full dose up to **Secondary IV Bag** - **Closed system safely handles hazardous fluids.**

Turn **Flow Selection Knob** to the West and pull full flush into the injector syringe

Turn **Flow Selection Knob** to the East & push flush up to **Secondary IV Bag** & close **Pink Roller Clamp**

Turn **Flow Selection Knob** to the North & open **Pink Roller Clamp** allowing IV Bag flush to drain into the chamber. Note: By controlling this flow rate, there can be a minimal wetting of the burette Chamber

Upon exhaustion of **Secondary IV Bag** fluid in the chamber, **AutoStart** mode commences to flush the line supporting vein patency

Drug Port Dose to Filled Chamber

Turn **Flow Selection Knob** to the South & add medication via injection through **Needless Drug Port Connector**

Turn **Flow Selection Knob** to the West and pull flush into injector syringe

Push syringe flush into the chamber. Check valve prevents flow into the Primary line. In the West position, Primary fluid may be allowed to fill the chamber to a prescribed level for dilution of the injected medication, or the **Flow Selection Knob** can be turned back to the South. **Note: The syringe may remain attached if hazardous material is in use for no exposure.**



Medication flow rate to the patient is managed by the White Roller Clamp establishing the prescribed drip rate in the Drip Chamber.

Flush of Dosed Secondary IV Bag

An injector syringe may be attached to the **Needless Drug Port Connector**. Turn the **Flow Control Knob** to the West and pull a flush into the injector

Turn **Flow Control Knob** to the East & push flush up to the **Secondary IV Bag**, then turn Knob to the North sending **Secondary IV Bag** flush to the chamber

Safety Concern	Risk	Design Feature
Air Embolism	Increased mortality	Weight of float shuttle valve ensures complete sealing of liquid to prevent air passage - operating to an extreme angle of 60° suitable for use in transport vehicles
Closed System	Infection / Contamination	Ready spike connections for Primary IV bag (typically 1,000 ml), Secondary IV bag, & a Drug Port for delivery of meds to the burette or to dose/flush the Secondary IV bag
Syringe Flush	Cross Contamination	Eliminates any need for a separate flushing syringe. Always ready flush pull from the AutoStart feed line
Secondary IV Bag Flush	Full Dose Delivery	A 50 ml Secondary Bag infusion may retain 7 ml of drug volume or 14% loss of the prescribed dose without adequate flush - SIS supports flushing Secondary IV bag
Catheter Flush	PIVC Failure	SIS AutoStart function automatically restarts IV flow after medication optimizing continuous flow conditions at the catheter preventing internal luminal occlusion
Medication Dilution	Interrupted Therapy	Burette filled to dilution specification from Primary IV bag, followed by precision dose at Drug Port or from a Secondary IV bag . After full dose delivery to burette, AutoStart initiates to accomplish full flush.
Protocol Adherence	Infection / Drug priority	Observational studies demonstrate a high incidence of non-conformance: 1) closed system to counter lack of hand washing 2) AutoStart ensures flushing protocol 3) Incorrect positioning of Secondary IV bag eliminated
Flushing Pressure	PIVC Failure	<p>2,000 mm³ pushed through a 20-gauge catheter (OD 1.1mm ID 0.85mm) in a vein (ID 2.5 to 4 mm) represents 3,500 mm of length or 11 feet</p> <p>2,000mm³ = 3,500mm x π x (0.85mm/2)² Volume of tube = length x π x radius²</p> <p>Best to flush through the burette for a steady drip rate</p> <p>Dispensing 2 ml during push-pause = 2,000 mm³ / push 5 times for a 10 ml dose</p> <p>3.5 meters 11 feet of liquid driven from the catheter tip per 2 ml push</p> <p>Bolus flush to catheter by syringe may vary based on size & push/pause protocol - AutoStart and gravity feed ensure low and constant flushing pressure</p>
Cybersecurity	Any electronic device	Relying on gravity-feed only, SIS eliminates all concerns and costs associated with electronic pumps
Documentation	Shift Change Handoff	Included sterile labels and marker ensure timely medication documentation; Establishes process ownership




35% to 50% of PIVCs fail, according to 2021 *Infusion Therapy Standards of Practice*, before their intended use is complete due to complications like phlebitis, infiltration, occlusion, or infection. Comparing infusion nurse-inserted catheters to generalist nurse insertions found lower rates of leakage, phlebitis, and infiltration (statistically significant), but still noted complications in 20-30% of cases even with skilled hands, (A 2005 study in the *Journal of Infusion Nursing* (Palefski & Stoddard))

How we are helping to improve outcomes :

Simplify IV fluid setups with a single, universal system that

- 1) Facilitates the correct static head for optimum gravity feed**
- 2) Eliminates the potential for air embolism in the IV line**
- 3) Adapted for multiple inputs all delivered with gentle drip flow**
- 4) Positions the drip chamber at eye level for caregiver awareness**
- 5) Re-establishes saline flow after any medicinal input without caregiver intervention for steady catheter flow to optimize vein patency**
- 6) Enables a ready minibag or syringe flush for complete drug delivery**
- 7) Emphasis on sterile labeling of each fluid event for quality of care**




Safer Infusion Systems solves your toughest infusion challenges – handling hazardous antineoplastic compounds. **SIS Oncology™**'s comprehensive materials approach takes caustic and light-sensitive fluids in stride, protecting the caregiver and patient while optimizing infusion success through gentle, gravity fluid delivery replacing electronic pumping to enhance vein patency. Note: When extreme precision is necessary, the SIS is entirely compatible with electronic pump devices.



Pigmentation of molded & extruded materials blocks UV & visible light

0.2-micron PTFE filter membrane stacked with carbon fiber membrane for vapor management




The Gravity Grand Slam

1 Hang single closed system SIS burette

Large saline + AutoStart feed 2 steady flush

3 Sterile input options



for IV infusion success:

Infection prevention

Promotion of vein preservation

Guarantee of proper static head gravity feed

Air embolism prevention with float valve

Ask for our white paper: *"Revitalizing Gravity-Fed IVs to Improve Vein Patency and Eliminate Caregiver Errors Through Standardization"*

Contact for distribution / licensing
(805) 377-1690
tcampbell@saferinfusion.com



Ver. 03-18-2025

Heart of AutoStart: Buoyant volume & float weight activate this shuttle valve for automatic gentle IV flow optimizing vein patency

